

FITC anti-mouse Ly-49C/F/H/I

Catalog No: #FC1042



Package Size: #FC1042-1 25Tests #FC1042-2 100Tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FITC anti-mouse Ly-49C/F/H/I
Brief Description	Syrian Hamster Monoclonal Antibody
Host Species	Syrian Hamster
Clonality	Monoclonal
Clone No.	s-14B11
Isotype	Syrian Hamster IgG1
Applications	FC
Species Reactivity	mouse
Conjugates	FITC
Target Name	Ly-49C/F/H/I
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only