

## PE anti-mouse CD3

Catalog No: #FC1096

Package Size: #FC1096-1 25Tests #FC1096-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	PE anti-mouse CD3
Brief Description	rat monoclonal antibody
Host Species	rat
Clonality	Monoclonal
Clone No.	s-17A2
Isotype	rat IgG2b, $\kappa$
Applications	FC
Species Reactivity	mouse
Conjugates	PE
Target Name	CD3
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Note: This product is for in vitro research use only