

Pacific Blue anti-mouse CD69

Catalog No: #FC1209



Package Size: #FC1209-1 25Tests #FC1209-2 100Tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Pacific Blue anti-mouse CD69 |
| Brief Description | Armenian Hamster Monoclonal Antibody |
| Host Species | Armenian Hamster |
| Clonality | Monoclonal |
| Clone No. | s-H1.2F3 |
| Isotype | Armenian Hamster IgG1 |
| Applications | FC |
| Species Reactivity | mouse |
| Conjugates | Pacific Blue |
| Target Name | CD69 |
| Formulation | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL. |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only