

Alexa Fluor 488 anti-mouse IL-4

Catalog No: #FC1223



Package Size: #FC1223-1 25Tests #FC1223-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Alexa Fluor 488 anti-mouse IL-4
Brief Description	Rat monoclonal antibody
Host Species	Rat
Clonality	Monoclonal
Clone No.	s-11B11
Isotype	Rat IgG1,k
Applications	FC
Species Reactivity	mouse
Conjugates	Alexa Fluor 488
Target Name	IL-4
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Note: This product is for in vitro research use only