

PE-Cy5 Mouse IgG1 κ Isotype Control

Catalog No: #FC1291



Package Size: #FC1291-1 25Tests #FC1291-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PE-Cy5 Mouse IgG1 κ Isotype Control
Clonality	Monoclonal
Clone No.	s-IB7.11
Applications	FC
Species Reactivity	Mouse
Conjugates	PE-Cy5
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only