

PE anti-human TNF- $\alpha$ 

Catalog No: #FC1303

Package Size: #FC1303-1 25Tests #FC1303-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	PE anti-human TNF- $\alpha$
Clonality	Monoclonal
Clone No.	s-F6
Applications	FC
Species Reactivity	human
Conjugates	PE
Target Name	TNF- $\alpha$
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 10  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Intracellular Cell Stimulation Cocktail (plus protein transport inhibitors)-stimulated peripheral blood mononuclear cells was fixed permeabilized and intracellularly stained with PE anti-human TNF- $\alpha$

Note: This product is for in vitro research use only