

## PE-Cy7 anti-human CD95

Catalog No: #FC1398

Package Size: #FC1398-1 25Tests #FC1398-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	PE-Cy7 anti-human CD95
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-DX2
Isotype	mouse IgG1, $\kappa$
Applications	FC
Species Reactivity	human
Conjugates	PE-Cy7
Target Name	CD95
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood lymphocytes stained with PE-Cy7 anti-human CD95 (colored histogram) or PE-Cy7 mouse IgG1  $\kappa$  (gray histogram)

Note: This product is for in vitro research use only