

PE anti-human/mouse CD44

Catalog No: #FC1403



Package Size: #FC1403-1 25 Tests #FC1403-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	PE anti-human/mouse CD44
Clonality	Monoclonal
Clone No.	s-IM7
Isotype	rat IgG2b, κ
Applications	FC
Species Reactivity	human/mouse
Conjugates	PE
Target Name	CD44
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with PE anti-human CD44 (colored histogram) or PE rat IgG2b κ gray histogram)

Note: This product is for in vitro research use only