

# Purified anti-human CD2

Catalog No: #FC1421



Package Size: #FC1421-1 100ug #FC1421-2 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                         |
|--------------------|-------------------------|
| Product Name       | Purified anti-human CD2 |
| Brief Description  | Rat Monoclonal Antibody |
| Host Species       | Rat                     |
| Clone No.          | s-LO-CD2a               |
| Species Reactivity | human                   |
| Storage            | Store at 4°C.           |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only