

Purified anti-human CD23

Catalog No: #FC1464



Package Size: #FC1464-1 100ug #FC1464-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human CD23
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clone No.	s-EBVCS-5
Species Reactivity	human
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only