

Purified anti-Mouse/Human CD44

Catalog No: #FC1484



Package Size: #FC1484-1 100ug #FC1484-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-Mouse/Human CD44
Brief Description	rat monoclonal antibody
Host Species	rat
Clone No.	s-IM7
Species Reactivity	Mouse/Human
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only