

Purified anti-human CD45RO

Catalog No: #FC1493



Package Size: #FC1493-1 100ug #FC1493-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|----------------------------|
| Product Name | Purified anti-human CD45RO |
| Brief Description | mouse monoclonal antibody |
| Host Species | mouse |
| Clone No. | s-UCHL-1 |
| Species Reactivity | human |
| Storage | Store at 4°C. |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only