

Purified anti-human CD152

Catalog No: #FC1533



Package Size: #FC1533-1 100ug #FC1533-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human CD152
Brief Description	human monoclonal antibody
Host Species	human
Clone No.	s-MDX-01-1000
Species Reactivity	human
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only