

Purified anti-human/mouse/rat CD278 (ICOS)

Catalog No: #FC1542



Package Size: #FC1542-1 100ug #FC1542-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human/mouse/rat CD278 (ICOS)
Brief Description	Armenian Hamster Monoclonal Antibody
Host Species	Armenian Hamster
Clone No.	s-C398.s-
Species Reactivity	human/mouse/rat
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only