

Purified anti-human HLA-A α B α C α HLA Class I α



Catalog No: #FC1550

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human HLA-A α B α C α HLA Class I α
Brief Description	mouse monoclonal antibody
Clone No.	s-W6/32
Species Reactivity	human
Storage	Store at 4°C.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only