

# Purified anti-human HLA-A $\alpha$ B $\alpha$ C $\alpha$ HLA Class I $\alpha$



Catalog No: #FC1550

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Purified anti-human HLA-A $\alpha$ B $\alpha$ C $\alpha$ HLA Class I $\alpha$
Brief Description	mouse monoclonal antibody
Clone No.	s-W6/32
Species Reactivity	human
Storage	Store at 4°C.

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only