

# Purified anti-human IL-12p40

Catalog No: #FC1558



Package Size: #FC1558-1 100ug #FC1558-2 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Purified anti-human IL-12p40
Brief Description	Human monoclonal antibody
Host Species	Human
Clone No.	s-ABT-874
Species Reactivity	human
Storage	Store at 4°C.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only