

Purified anti-human IL-12p70

Catalog No: #FC1559



Package Size: #FC1559-1 100ug #FC1559-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human IL-12p70
Brief Description	rat monoclonal antibody
Host Species	rat
Clone No.	s-20-c2
Species Reactivity	human
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only