

Purified anti-human iNKT cell (TCR Vα24-Jα18)

Catalog No: #FC1564



Package Size: #FC1564-1 100ug #FC1564-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human iNKT cell (TCR Vα24-Jα18)
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clone No.	s-6B11
Species Reactivity	human
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only