

Purified anti-mouse CD8b.2

Catalog No: #FC1585



Package Size: #FC1585-1 100ug #FC1585-2 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Purified anti-mouse CD8b.2
Brief Description	rat monoclonal antibody
Host Species	rat
Clone No.	s-53-5.8
Species Reactivity	mouse
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only