

# Purified anti-Mouse/Human CD49d

Catalog No: #FC1599



Package Size: #FC1599-1 100ug #FC1599-2 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Purified anti-Mouse/Human CD49d
Brief Description	rat monoclonal antibody
Host Species	rat
Clone No.	s-PS/2
Species Reactivity	Mouse/Human
Storage	Store at 4°C.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only