

Purified anti-mouse CD183

Catalog No: #FC1632



Package Size: #FC1632-1 100ug #FC1632-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse CD183
Brief Description	Armenian Hamster Monoclonal Antibody
Host Species	Armenian Hamster
Clone No.	s-CXCR3-173
Species Reactivity	mouse
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only