

Purified anti-mouse H-2Kb

Catalog No: #FC1657



Package Size: #FC1657-1 100ug #FC1657-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse H-2Kb
Brief Description	Mouse monoclonal antibody
Host Species	Mouse
Clone No.	s-AF6-88.5
Species Reactivity	mouse
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only