

Purified anti-mouse H-2Kd (MHC Class I)

Catalog No: #FC1658



Package Size: #FC1658-1 100ug #FC1658-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Purified anti-mouse H-2Kd (MHC Class I) |
| Brief Description | mouse monoclonal antibody |
| Host Species | mouse |
| Clone No. | s-A4C8.1-Do9 |
| Species Reactivity | mouse |
| Storage | Store at 4°C. |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only