

Purified anti-mouse H-2Kd (MHC Class I)

Catalog No: #FC1658



Package Size: #FC1658-1 100ug #FC1658-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse H-2Kd (MHC Class I)
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clone No.	s-A4C8.1-Do9
Species Reactivity	mouse
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only