

# Purified anti-mouse IFN $\gamma$

Catalog No: #FC1662



Package Size: #FC1662-1 100ug #FC1662-2 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Purified anti-mouse IFN $\gamma$
Brief Description	rat monoclonal antibody
Host Species	rat
Clone No.	s-XMG1.2
Species Reactivity	mouse
Storage	Store at 4°C.

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only