

Purified anti-mouse Ly-6C

Catalog No: #FC1668



Package Size: #FC1668-1 100ug #FC1668-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse Ly-6C
Brief Description	Armenian Hamster Monoclonal Antibody
Host Species	Armenian Hamster
Clone No.	s-7B10
Species Reactivity	mouse
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only