

Purified anti-mouse Ly-6G/Ly-6C

Catalog No: #FC1670



Package Size: #FC1670-1 100ug #FC1670-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse Ly-6G/Ly-6C
Brief Description	rat monoclonal antibody
Host Species	rat
Clone No.	s-RB6-8C5
Species Reactivity	mouse
Storage	Store at 4°C.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only