

Purified anti-mouse TCR $\gamma\delta$

Catalog No: #FC1674



Package Size: #FC1674-1 100ug #FC1674-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-mouse TCR $\gamma\delta$
Brief Description	Hamster Monoclonal Antibody
Host Species	Hamster
Clone No.	s-UC7-13D5
Species Reactivity	mouse
Storage	Store at 4°C.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only