

Purified anti-human/mouse TGF- β 1 2 3

Catalog No: #FC1676



Package Size: #FC1676-1 100ug #FC1676-2 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|----------------------------------------------|
| Product Name | Purified anti-human/mouse TGF- β 1 2 3 |
| Brief Description | mouse monoclonal antibody |
| Host Species | mouse |
| Clone No. | s-1D11.16.8 |
| Species Reactivity | human/mouse |
| Storage | Store at 4°C. |

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only