

# Purified anti-human/mouse TGF- $\beta$ 1 2 3

Catalog No: #FC1676



Package Size: #FC1676-1 100ug #FC1676-2 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Purified anti-human/mouse TGF- $\beta$ 1 2 3
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clone No.	s-1D11.16.8
Species Reactivity	human/mouse
Storage	Store at 4°C.

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only