

# Purified anti-mouse V $\gamma$ 2 TCR

Catalog No: #FC1679



Package Size: #FC1679-1 100ug #FC1679-2 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Purified anti-mouse V $\gamma$ 2 TCR
Brief Description	Armenian Hamster Monoclonal Antibody
Host Species	Armenian Hamster
Clone No.	s-UC3-10A6
Species Reactivity	mouse
Storage	Store at 4°C.

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only