

VPRBP rabbit pAb

Catalog No: #62018

Package Size: #62018-1 50ul #62018-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	VPRBP rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of VPRBP at Human/Mouse
Immunogen Description	Synthesized peptide derived from human VPRBP
Target Name	VPRBP DCAF1 KIAA0800 RIP
Uniprot	Q9Y4B6
GeneID	9730
Calculated MW	166kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1αO 500-2000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only