

# Bovine Growth Hormone (GH) ELISA Kit

Catalog No: #EK20775



Package Size: #EK20775-1 48T #EK20775-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine Growth Hormone (GH) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Bovine

## Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:28pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only