

# Bovine Immunoglobulin A (IgA) ELISA Kit

Catalog No: #EK20787



Package Size: #EK20787-1 48T #EK20787-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine Immunoglobulin A (IgA) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Bovine

## Application Details

Detect Range: 31.2-2000  $\mu$ g/mL

Sensitivity: 12.6  $\mu$ g/mL

Sample Type: Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50  $\mu$ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only