

LMOD2 rabbit pAb

Catalog No: #62019

Package Size: #62019-1 50ul #62019-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LMOD2 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA;IHC
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of LMOD2 at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human LMOD2
Target Name	LMOD2
Uniprot	Q6P5Q4
GeneID	442721
Calculated MW	60kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-2000;IHC 1:50-300; ELISA 2000-20000

Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung.

Note: This product is for in vitro research use only