

# Human Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) ELISA Kit



Catalog No: #EK20777

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20777-1 48T #EK20777-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.312-20ng/mL
Sensitivity:0.125ng/mL
Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only