

SLC38A5 rabbit pAb

Catalog No: #62020



Package Size: #62020-1 50ul #62020-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC38A5 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of SLC38A5 at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human SLC38A5
Target Name	SLC38A5 JM24 SN2 SNAT5 PP7194
Uniprot	Q8WUX1
GeneID	92745
Calculated MW	52kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only