

IRG1 Rabbit Monoclonal Antibody

Catalog No: #62023



Package Size: #62023-1 50ul #62023-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IRG1 Rabbit Monoclonal Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity Purification
Applications	WB,IP
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human IRG1
Other Names	CAD; ACOD1
Uniprot	A6NK06,P54987
GeneID	73024916365
Calculated MW	53kDa
Formulation	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000, IP 1:20

Note: This product is for in vitro research use only