

## Ly6g Rabbit Monoclonal Antibody

Catalog No: #62028



Package Size: #62028-1 50ul #62028-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ly6g Rabbit Monoclonal Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity Purification
Applications	IHC-P,FCM,IF-P
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of mouse Ly6g
Other Names	Gr1; Gr-1; Ly-6G
Uniprot	P35461
GeneID	546644
Calculated MW	14kDa
Formulation	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Storage	Store at -20°C/1 year

## Application Details

IHC-P/IF-P 1:500-1:2000, FCM 1:500

Note: This product is for in vitro research use only