

NSUN5 Rabbit pAb

Catalog No: #62026



Package Size: #62026-1 50ul #62026-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NSUN5 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC/IF
Species Reactivity	Human
Immunogen Description	KLH conjugated synthetic peptide derived from human NSUN5
Other Names	FLJ10267; MGC986; NOL1; NOL1-related protein; NOL1/NOP2/Sun domain family member 5; NOL1R; NOP2/Sun domain family, member 5; NOP2/Sun domain family, member 5A; NSUN5; NSUN5_HUMAN; p120; Putative methyltransferase NSUN5; WBSCR20; WBSCR20A; Williams Beuren syndrome chromosome region 20A; Williams-Beuren syndrome chromosomal region 20A protein; Williams-Beuren syndrome critical region protein 20 copy A.
Uniprot	Q96P11
GeneID	55695
Calculated MW	47 kDa
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, ICC/IF=1:100-500

Note: This product is for in vitro research use only