

# Recombinant Myc-Tag Monoclonal Antibody(2G4), HRP Conjugated

Catalog No: #T401

Package Size: #T401-1 50ul #T401-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Myc-Tag Monoclonal Antibody(2G4), HRP Conjugated
Purification	affinity purified by Protein A
Applications	WB,ELISA
Species Reactivity	Species independent
Conjugates	HRP
Formulation	1Γ PBS (pH7.4) oO 0.05% proclin300oO 50% Glycerol
Storage	Store at -20°C/1 year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:30000-1:100000, ELISA 1:5000-1:10000

Note: This product is for in vitro research use only