

# Recombinant Flag-Tag Monoclonal Antibody(3O2), HRP Conjugated



Catalog No: #T402

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #T402-1 50ul #T402-2 100ul

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Flag-Tag Monoclonal Antibody(3O2), HRP Conjugated
Purification	affinity purified by Protein A
Applications	WB,ELISA
Species Reactivity	Species independent
Conjugates	HRP
Formulation	1Γ PBS (pH7.4) oO 0.05% proclin300oO 50% Glycerol
Storage	Store at -20°C/1 year. Avoid repeated freeze/thaw cycles.

## Application Details

WB 1:30000-1:100000, ELISA 1:5000-1:10000

Note: This product is for in vitro research use only