

Recombinant HA-Tag Monoclonal Antibody(12H2), HRP Conjugated



Catalog No: #T403

Orders: order@signalwayantibody.com

Package Size: #T403-1 50ul #T403-2 100ul

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Recombinant HA-Tag Monoclonal Antibody(12H2), HRP Conjugated |
| Purification | affinity purified by Protein A |
| Applications | WB,ELISA |
| Species Reactivity | Species independent |
| Conjugates | HRP |
| Formulation | 1Γ PBS (pH7.4) oO 0.05% proclin300oO 50% Glycerol |
| Storage | Store at -20°C/1 year. Avoid repeated freeze/thaw cycles. |

Application Details

WB 1:30000-1:100000, ELISA 1:5000-1:10000

Note: This product is for in vitro research use only