

# Recombinant His-Tag Monoclonal Antibody(14W7), HRP Conjugated

Catalog No: #T404

Package Size: #T404-1 50ul #T404-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Recombinant His-Tag Monoclonal Antibody(14W7), HRP Conjugated |
| Purification       | affinity purified by Protein A                                |
| Applications       | WB,ELISA  |
| Species Reactivity | Species independent   |
| Conjugates         | HRP   |
| Formulation        | 1Γ PBS (pH7.4) oO 0.05% proclin300oO 50% Glycerol             |
| Storage            | Store at -20°C/1 year. Avoid repeated freeze/thaw cycles.     |

## Application Details

WB 1:30000-1:100000, ELISA 1:5000-1:10000

Note: This product is for in vitro research use only