

## EphA2 (Phospho Ser897) rabbit pAb

Catalog No: #62027



Package Size: #62027-1 50ul #62027-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	EphA2 (Phospho Ser897) rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of Human Mouse EphA2 (phospho-Ser897)
Immunogen Description	Synthesized phospho peptide around human EphA2 (Ser897)
Target Name	EPHA2 ECK
Uniprot	P29317
GeneID	1969
Calculated MW	105kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only