

Vinculin mouse mAb

Catalog No: #62029



Package Size: #62029-1 50ul #62029-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vinculin mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Purification	Protein G
Applications	WB;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of Vinculin protein.
Immunogen Description	Synthesized peptide derived from human Vinculin.
Target Name	VCL
Other Names	VCL;Vinculin;Metavinculin
Uniprot	P18206
GeneID	7414
Calculated MW	124kD
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000. IF 1:100-500. ELISA 1:1000-5000

Note: This product is for in vitro research use only