

## Anti-PINK1 (Phospho-S228) Antibody

Catalog No: #62030



Package Size: #62030-1 50ul #62030-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | Anti-PINK1 (Phospho-S228) Antibody   |
| Host Species          | Rabbit   |
| Clonality             | Polyclonal   |
| Purification          | The antibody was purified by immunogen affinity chromatography.  |
| Applications          | WB, IH   |
| Species Reactivity    | H, M, R  |
| Immunogen Description | KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S228 of human PINK1 protein. |
| Other Names           | Serine/threonine-protein kinase PINK1 mitochondrial; BRPK; PTEN-induced putative kinase protein 1          |
| Uniprot               | Q9BXM7   |
| GeneID                | 65018  |
| Formulation           | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.  |
| Storage               | store at -20°C for one year  |

## Application Details

WB (1/500 - 1/1000), IH (1/50 - 1/100)

Note: This product is for in vitro research use only