

EOMES Polyclonal Antibody

Catalog No: #62031

Package Size: #62031-1 50ul #62031-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EOMES Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse
Specificity	EOMES Polyclonal Antibody detects endogenous levels of EOMES protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human EOMES.
Target Name	EOMES
Other Names	EOMES;TBR2;Eomesodermin homolog;T-box brain protein 2;T-brain-2;TBR-2
Uniprot	O95936
GeneID	8320
Calculated MW	73kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000. ELISA: 1:10000.

Images



Western blot analysis of lysate from HT29 cells, using EOMES antibody.

Note: This product is for in vitro research use only