

Mouse Cytokeratin Fragment Antigen 21-1 (CYFRA21-1) ELISA Kit

Catalog No: #EK20786

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20786-1 48T #EK20786-2 96T

Description

Product Name	Mouse Cytokeratin Fragment Antigen 21-1 (CYFRA21-1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:6.38pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only